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- Gross, A.S. (2001). Best practice in therapeutic drug monitoring. *J Clin Pharmacol* 52, 5S–10S.
- American Society of Health-System Pharmacists. (2008). Pharmacogenetics, therapeutic drug management offer new insights. *Pharmacy News* 1 September 2008. Available at: <http://www.ashp.org/menu/News/PharmacyNews/NewsArticle.aspx?id=2902>
- George, C.H., Barberini-Jammaers, S.R., Muller, C.T. (2008). Refocusing therapeutic strategies for cardiac arrhythmias: defining viable molecular targets to restore cardiac ion flux. *Expert Opin Ther Patents* 18, 1–19.
- Campbell, T.J., Williams, K.M. (1998). Therapeutic drug monitoring: antiarrhythmic drugs. *Br J Clin Pharmacol* 46, 307–19.
- Digoxin: serious drug interactions. (2010). *Prescribe Int* 19, 68–70.
- Heard, K. (2004). Therapeutic drugs. In: *Medical Toxicology*, 3rd edition. Dart, R.C. ed. Lippincott Williams and Wilkins, Philadelphia, PA.
- Connolly, S.J., Kates, R.E. (1982). Clinical pharmacokinetics of N-acetylprocainamide *Clin Pharmacokinet* 7, 206–20.
- Smellie, W.S., Coleman, J.J. (2007). Pitfalls of testing and summary of guidance on safety monitoring with amiodarone and digoxin. *BMJ* 334, 312–5.
- Cobo, A., Martín-Suarez, A., Calvo, M.V., Domínguez-Gil, A., de Gatta, M.M. (2010). Clinical repercussions of analytical interferences due to aldosterone antagonists in digoxin immunoassays: an assessment. *Ther Drug Monit* 32, 169–76.
- Dasgupta, A. (2006). Therapeutic drug monitoring of digoxin: impact of endogenous and exogenous digoxin-like immunoreactive substances. *Toxicol Rev* 25, 273–81.
- Dasgupta, A. (2008). Herbal supplements and therapeutic drug monitoring: focus on digoxin immunoassays and interactions with St. John's wort. *Ther Drug Monit* 30, 212–7.
- WHO (2009). Epilepsy. Fact sheet No. 999. Available at: <http://www.who.int/mediacentre/factsheets/fs999/en/>
- Tomson, T., Dahl, M.L., Kimland, E. (2007). Therapeutic monitoring of antiepileptic drugs for epilepsy. *Cochrane Database Syst Rev* CD002216.
- Patsalos, P.N., Berry, D.J., Bourgeois, B.F., Cloyd, J.C., Glauser, T.A., Johannessen, S.I., Leppik, I.E., Tomson, T., Perucca, E. (2008). Antiepileptic drugs--best practice guidelines for therapeutic drug monitoring: a position paper by the subcommission on therapeutic drug monitoring, ILAE Commission on Therapeutic Strategies. *Epilepsia* 49, 1239–76.
- Tomson, T., Johannessen, S.I. (2000). Therapeutic monitoring of the new antiepileptic drugs. *Eur J Clin Pharmacol* 55, 697–705.
- Vajda, F.J. (2007). Monitoring antiepileptic drug therapy with serum level measurements. *Med J Aust* 187, 581.
- Striano, S., Striano, P., Capone, D., Pisani, F. (2008). Limited place for plasma monitoring of new antiepileptic drugs in clinical practice. *Med Sci Monit* 14, RA173–8.
- WHO (2011). Asthma facts. Available at: <http://www.who.int/respiratory/asthma/scope/en/index.html>
- Lippi, G., Salvagno, G.L., Brocco, G., Guidi, G.C. (2007). Appropriateness of theophylline plasma levels. *Arch Pathol Lab Med* 131, 1166.
- Raebel, M.A., Carroll, N.M., Andrade, S.E., Chester, E.A., Lafata, J.E., Feldstein, A., Gunter, M.J., Nelson, W.W., Simon, S.R., Chan, K.A., Davis, R.L., Platt, R. (2006). Monitoring of drugs with a narrow therapeutic range in ambulatory care. *Am J Manag Care* 12, 268–74.
- Langen, U., Knopf, H., Melchert, H.U. (2006). Drug therapy of asthma. Results of the 1998 German National Health Interview and Examination Survey. *Bundesgesundheitsblatt Gesundheitsforschung Gesundheitsschutz* 49, 903–10.
- Merikangas, K.R., Jin, R., He, J.-P., Kessler, R.C., Lee, S., Sampson, N.A., Viana, M.C., Andrade, L.H., Hu, C., Karam, E.G., Ladea, M., Medina-Mora, M.E., Ono, Y., Posada-Villa, J., Sagor, R., Wells, J.E., Zarkov, Z. (2011). Prevalence and correlates of bipolar spectrum disorder in the world mental health survey initiative. *Arch Gen Psychiatry* 68, 241–51.
- Bowden, C.L. (2009). Anticonvulsants in bipolar disorders: current research and practice and future directions. *Bipolar Disord* 11 Suppl 2, 20–33.
- Center for Quality Assessment and Improvement in Mental Health. Available at: http://www.cqaimh.org/measure_LS.html
- Bentué-Ferrer, D., Tribut, O., Verdier, M.C., Suivi Thérapeutique Pharmacologique de la Société Française de Pharmacologie et de Thérapeutique. (2010). Therapeutic drug monitoring of valproate. *Therapie* 65, 233–40.
- Keck, P.E. Jr, McIntyre, R.S., Shelton, R.C. (2007). Bipolar depression: best practices for the outpatient. *CNS Spectr* 12 Suppl 20, 1–14.
- Paton, C., Barnes, T.R., Shingleton-Smith, A., McAllister-Williams, R.H., Kirkbride, J., Jones, P.B., McIntyre, S.; POMH-UK project team. (2010). Lithium in bipolar and other affective disorders: prescribing practice in the UK. *J Psychopharmacol* 24, 1739–46.
- Colom, F., Vieta, E.M., Martínez-Arán, A., Reinares, M., Benabarre, A., Gastó, C. (2000). Clinical factors associated with treatment noncompliance in euthymic bipolar patients. *J Clin Psychiatry* 61, 549–55.
- Rybak, M.J., Lomaestro, B.M., Rotschafer, J.C., Moellering, R.C. Jr, Craig, W.A., Billeter, M., Daloviso, J.R., Levine, D.P. (2009). Therapeutic monitoring of vancomycin in adults: summary of consensus recommendations from the American Society of Health-System Pharmacists, the Infectious Diseases Society of America, and the Society of Infectious Diseases Pharmacists. *Pharmacotherapy* 29, 1275–9.
- Global RPh. Therapeutic drug levels. Available at: <http://www.globalrph.com/levels.htm>
- Jelassi, M.L., Benmoudou, A., Lefeuve, S., Mainardi, J.L., Billaud, E.M. (2011). Level of evidence for therapeutic drug monitoring of vancomycin. *Therapie* 66, 29–37.
- Lopez-Novoa, K.M., Quiros, Y., Vicente, L., Moreales, A.I., Lopez-Hernandez, F.J. (2011). New insights into the mechanism of aminoglycoside nephrotoxicity: an integrative point of view. *Kidney Int* 79, 33–45.
- Conil, J.M., Georges, B., Breden, A., Segonds, C., Lavit, M., Seguin, T., Coley, N., Samii, K., Chabanon, G., Houin, G., Saivin, S. (2006). Increased amikacin dosage requirements in burn patients receiving a once-daily regimen. *Int J Antimicrob Agents* 28, 226–30.
- Dolton, M., Xu, H., Cheong, E., Maitz, P., Kennedy, P., Gottlieb, T., Buono, E., McLachlan, A.J. (2010). Vancomycin pharmacokinetics in patients with severe burn injuries. *Burns* 36, 469–76.
- Le Floch, R., Arnould, J.F., Pilonget, A., Dally, E., Naux, E. (2010). Antimicrobial blood concentrations in burns. A five years' retrospective survey. *Pathol Biol (Paris)* 58, 137–43.
- Taccone, F.S., Laterre, P.F., Spapen, H., Dugernier, T., Delattre, I., Layeux, B., De Backer, D., Wittebole, X., Wallemacq, P., Vincent, J.L., Jacobs, F. (2010). Revisiting the loading dose of amikacin for patients with severe sepsis and septic shock. *Crit Care* 14, R53.
- del Mar Fernández de Gatta García, M., Revilla, N., Calvo, M.V., Domínguez-Gil, A., Sánchez Navarro, A. (2007). Pharmacokinetic/pharmacodynamic analysis of vancomycin in ICU patients. *Intensive Care Med* 33, 279–85.
- Rea, R.S., Capitano, B., Bies, R., Bigos, K.L., Smith, R., Lee, H. (2008). Suboptimal aminoglycoside dosing in critically ill patients. *Ther Drug Monit* 30, 674–81.
- Budde, K., Glande, P. (2011). Pharmacokinetic principles of immunosuppressive drugs. *Clinical Pharmacology and Therapeutics* 90, 316–22.
- Kuypers, D.R. et al. (2008). Influence of interactions between immunosuppressive drugs on therapeutic drug monitoring. *Annals of Transplantation* 13, 11–18.
- Levitzy, J. et al. (2011) Next level of immunosuppression: drug/immune monitoring. *Liv Trans* 17, 60–65.
- Schiff, J., Cole, E., Cantaorvich, M. (2007) Therapeutic monitoring of calcineurin inhibitors for the nephrologist. *Liv Clinical Journal of the American Society of Nephrology* 2, 374–384.
- Levy, G. et al. (2002) Improved clinical outcome for liver transplant recipients using cyclosporine monitoring based on 2-hr post-dose levels (C2). *Transplantation* 73, 953–959.
- Wallemacq, P. et al. (2009). Opportunities to optimize tacrolimus therapy in solid organ transplantation. *Therapeutic drug monitoring* 31, 139–152.
- Kuypers, D. et al. (2012). Consensus report on therapeutic drug monitoring of mycophenolic acid in solid organ transplantation. *Clinical Journal of the American Society of Nephrology* 5, 341–358.
- Holt, D.W. et al. (2003). Therapeutic monitoring of sirolimus: its contribution to optimal prescription. *Transplantation Proceedings* 35, 157–161.
- Rial, M. et al. (2010). Individualizing early use of sirolimus in renal transplantation. *Transplantation Proceedings* 42, 4518–4525.
- Sumethkul, V. et al. (2010). Long-term outcome of very early cyclosporine minimization and de novo everolimus therapy in kidney transplant recipients: a pharmacokinetic approach. *Transplantation Proceedings* 42, 4040–4042.
- Casanovas, T., Argudo, A., Pena-Cala, M.C. (2011). Everolimus in clinical practice in long-term liver transplantation: an observational study. *Transplantation Proceedings* 43, 2216–2219.
- van Rossum, H.H., de Fijter, J.W., van Pelt, J. (2010). Pharmacodynamic monitoring of calcineurin inhibition therapy: principles, performance and perspectives. target ranges of mycophenolic acid exposure and drug-related adverse events: a 5-year, open-label, prospective, clinical follow-up study in renal allograft recipients. *Therapeutic Drug Monitoring* 32, 3–10.

51. Jorga, A., Holt, D.W., Johnson, A. (2004). Therapeutic drug monitoring of cyclosporine. *Transplantation Proceedings* 36, 396–403.
52. Ekberg, H. et al. (2010). Cyclosporine, Tacrolimus and sirolimus retain their distinct toxicity profiles despite low doses in the Symphony Study. *Nephrology, Dialysis, Transplantation* 25, 2004–2010.
53. Ciancio, G. et al. (2006). A randomized long-term trial of tacrolimus/sirolimus versus tacrolimus/mycophenolate versus cyclosporine/sirolimus in renal transplantations: three-year analysis. *Transplantation* 81, 845–852.
54. Gaston, R.S. et al. (2009). Fixed or controlled-dose mycophenolate mofetil with standard or reduced dose calcineurin inhibitors. *American Journal of Transplantation* 9, 1607–1619.
55. van Gelder, T. et al. (2010). Renal transplant patients at high risk of acute rejection benefit from adequate exposure to mycophenolic acid. *Transplantation* 89, 595–599.
56. Musuamba, F.T., Rousseau, A., Bosmans, J.L., Senessaël, J.J., Cumps, J., Marquet, P., Wallemacq, P., Verbeeck, R.K. (2009). Limited sampling models and Bayesian estimation for mycophenolic acid under the curve prediction in stable renal transplant patients co-medicated with cyclosporin or sirolimus. *Clinical Pharmacokinetics* 48, 745–758.
57. Wu, M.J., et al. (2008). Impact of variability of sirolimus trough levels on chronic allograft nephropathy. *Transplantation Proceedings* 40, 2202–2205.
58. Lebranchua, Y., Snaouidjb, R., Toupancec, O., Weesteld, P.F., de Lignye, B.H., Buchlera, M., Rerollef, A., Thierryg, B., Moulinh, B., Subrai, J.-F., Deteixj, P., Pogampk, P.L., Finzil, L., Etiennem, I. (2012). Five-Year results of a randomized trial comparing de novo sirolimus and cyclosporine in renal transplantation. *American Journal of Transplantation*.
59. Kovarik, J.M. et al. (2001) Longitudinal assessment of everolimus in de novo renal transplant recipients over the first post-transplant year: pharmacokinetics, exposure-response relationships, and influence of cyclosporine. *Clinical Pharmacology and Therapeutics* 69, 48–56.
60. White, S., Wong, S.H.Y. (1998). Standards of laboratory practice: analgesic drug monitoring. *Clin Chem* 44, 1110–23.
61. Niscola, P., Scaramucci, L., Vischini, G., Giovannini, M., Ferrannini, M., Massa, P., Tatangelo, P., Galletti, M., Palumbo, R. (2010). The use of major analgesics in patients with renal dysfunction. *Curr Drug Targets* 11, 752–8.
62. American Society of Anesthesiologists Task Force on Chronic Pain Management and the American Society of Regional Anesthesia and pain Medicine. (2010). Practice guidelines for chronic pain management. An updated report by the American Society of Anesthesiologists Task Force on Chronic Pain Management and the American Society of Regional Anesthesia and pain Medicine. *Anesthesiology* 112, 810–33.
63. Lee, T.J. (2010). Pharmacologic treatment for low back pain: one component of pain care. *Phys Med Rehabil Clin N Am* 21, 793–800.

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© 2012 Roche Diagnostics.
Roche Diagnostics Deutschland GmbH
Roche Diagnostics (Schweiz) AG
Roche Diagnostics GmbH